

# CLEANING EQUIPMENT

## FOUR HEAD CLEANING CAN

### 18 LITER (4.75 GAL) STAINLESS STEEL CLEANING CAN

Create your own custom cleaning can with the valves you use. Stainless steel containers should never be used with cleaning chemicals containing chlorine. Always verify ingredients before use.



201-420



201-420 can top view

PART NO.	DESCRIPTION
201-420	Complete with any four valves installed, specify valve part numbers




729-005

**D SYSTEM**  
Anheuser-Busch  
Beck's  
Coors Miller  
Yuengling



729-001

**S SYSTEM**  
Heineken  
Scottish and  
Newcastle



729-003

**A SYSTEM**  
Warsteiner  
Paulaner  
Spaten



729-004

**G SYSTEM**  
Boddington



729-006

**M SYSTEM**  
Schneider  
Aventinus  
Einbecker



729-002

**U SYSTEM**  
Guinness  
Harp

## DIRECT DRAW CLEANER

This simple cleaner and sterilizer is ideal for cleaning direct draw systems. Includes a shatter-proof plastic bottle, HP-300 pump and faucet cleaning attachment. Complete kit includes a faucet brush, faucet wrench, chemical and instructions.

PART NO.	DESCRIPTION
CK-1001	1 qt - cleaner only
CK-1002	2 qt - cleaner only
CK-1100	1 qt - complete kit
CK-1200	2 qt - complete kit
M5-801147-CK	5 L - complete kit

## HYGIENE METER

The hand-held ATP Hygiene Meter measures the amount of bacteria on a surface. It only takes 15 seconds for a line cleaner to use this device to measure beer lines, faucets, etc. before and after cleaning.

PART NO.	DESCRIPTION
ATP-100	ATP hygiene meter
ATP-SWAB	Pack of 100 swabs (NOTE: Must be refrigerated)

## PH TESTING

Safety first! Always test rinse water to ensure a pH of 7- 8.

PART NO.	DESCRIPTION
PH-METER	Digital pH meter with water temperature compensation
PH-PAPER	50' roll pH test paper



M5-801147-CK



CK-1100



ATP-SWAB



ATP-100



PH-METER



PH-PAPER