


# Adaptor Plate Cutout Dimensions












19194



For use only with items 19192, 19188, and 19198 adaptor plates

19198



ITEM #	OPENINGS	CASE LOT
19190	With two 6 3/8" (16.2 cm) inset holes – fits 78164 (4½ qt/3.9 L)	4
19192	With two 8 3/8" (21.3 cm) inset holes – fits 78184 (7¼ qt/6.9 L)	4
19194	Adaptor ring fits 19192, 19188, and 19198 adaptor plates – fits 78184 (7¼ qt/6.9 L)	2
72227	With one 8 3/8" (21.3 cm) and one 10 3/8" (26.4 cm) inset holes – fits 78184 (7¼ qt/6.9 L) & 78204 (11 qt/10.4 L)	1
19193	With one 4 7/8" (12.4 cm) and two 6 3/8" (16.2 cm) inset holes – fits 78154 (2½ qt/2.4 L) & two 78164 (4½ qt/3.9 L)	4
72228	With three 6 1/2" (16.5 cm) inset holes – fits 78164 (4½ qt/3.9 L) – fits thermoset wells only	1
19188*	With three 8 3/8" (21.3 cm) inset holes – fits 78184 (7¼ qt/6.9 L) – for use with items 72050/72051	1
19195	With six 4 1/4" (10.8 cm) holes – fits 78710 (1¼ qt/1.2 L)	4
19189	With one 10 1/2" (26.8 cm) hole – fits 78204 (11 qt/10.4 L)	4
19198	With one 8 3/8" (21.3 cm) hole – fits 78184 (7¼ qt/6.9 L)	2

\*Fits ½ rethermalizers and ¼ drop-ins only